

# Lift Aid Equipment Determination Grid

(Recommended example)

Resident Dependency	Height Adjustable	Full Sling Lift	Stand Assist Lift	Lift Walker	Stand Assist Aid	Gait Belt with Handles	Friction Reducing
0 - Independent	Recommended	Never	Rarely	Rarely	Rarely	Rarely	Rarely
1 – Supervision	Recommended	Rarely	Occasionally	Occasionally	Normally	Occasionally	Rarely
2 – Limited Assistance	Strongly Recommended	Rarely	Normally	Normally	Normally	Normally	Normally
3 – Extensive Assistance	Required	Normally	Normally	Normally	Rarely	Rarely	Always
4 – Total Dependence	Required	Always	Never	Never	Never	Never	Always

For a typical resident with the dependency status classification as shown, this grid indicates normal equipment requirements to conduct a safe transfer. Some residents